

DMB COPY

DETROIT-MACOMB HOSPITAL CORPORATION

30-02005

7815 EAST JEFFERSON / DETROIT, MICHIGAN 48214 / TELEPHONE 821-6000

June 10, 1985

Mr. D. J. Sreniawski, Chief
Nuclear Materials Safety
Section 2
U.S. Nuclear Regulatory Commission
Region III
799 Roosevelt Road
Glen Ellyn, IL 60137

Ref.: License N:
21-01190-05

Dear Mr. Sreniawski:

In response to your letter of May 17, 1985, I submit the following information to correct the non-compliance items indicated in your letter:

1. The Calcheck Kit will be used in accordance with the procedures established by the manufacturer, as specified in License Condition N:16.
2. The accuracy test for the dose calibrator will be performed annually as specified in License Condition N:17.

To assure compliance to these items, the Radiation Safety Officer will personally supervise the test performances until operation personnel are familiar with how the tests are to be run and the expected results. The first of these tests were satisfactorily performed on Friday, May 31, 1985. The results of these tests are enclosed.

If additional information is required, please contact Mr. C. J. Damico, Radiation Safety Officer for the Corporation at (313) 573-5237.

8507150442 850610
REG3 LIC30
21-01190-05 PDR

/sn

Sincerely,

Jack Ryan, M.D.
President

RECEIVED
JUL 01 1985
REGION III

Attachments

cc: Messrs. C. J. Damico, Gareth Mitchell, Richard T. Young
D. L. Otto, H.D.

IE07

JUL 1 1985

DETROIT MEMORIAL HOSPITAL

1420 ST. ANTOINE / DETROIT, MICHIGAN 48226 / 225-5000

SOUTH MACOMB HOSPITAL

11800 TWELVE MILE / WARREN, MICHIGAN 48093 / 573-5000

DATA SHEET #1: (To be completed only ONCE)*

Kit Calibration

TESTS SUPERVISED BY C. J. DAMICO, R.S.O.

All readings must be taken at lowest range setting available and converted to mCi units.

Performed
5/31/85

TUBES A	DISPLAYED ACTIVITY B	CALIBRATION FACTORS C
Black Only	107 mCi	1.00
Black Only	107 mCi	
Black Only	107 mCi	1.86
Black & Red	57.4 mCi	
Black Only	107 mCi	3.43
Black & Orange	31.2 mCi	
Black Only	107 mCi	13.67
Black & Yellow	7.83 mCi	
Black Only	107 mCi	49.31
Black & Green	2.17 mCi	
Black Only	107 mCi	139.0
Black & Blue	0.77 mCi	
Black Only	107 mCi	449.6
Black & Purple	0.238 mCi	

SOURCE CONFIGURATION

Syringe

Vial

*Or following repair of dose calibrator or Cellcheck Kit. In all instances these factors can only be determined following proof of activity linearity by standard techniques. KEEP THIS FORM FOR FUTURE REFERENCE!

DATA SHEET #2 (to be completed each quarter)

Dose Calibrator Activity Linearity Check

Dose Calibrator Redx Date 5/31/85
Model Assayer 1 Technologist F. Allen
Source Configuration Vial (must be same as on Data Sheet #1)

All readings must be taken at lowest range setting available and converted to mCi units.

A	B	C	D
TUBE COLOR	DISPLAYED ACTIVITY	CALIBRATION FACTOR	PRODUCT OF B X C
Black Only:	107 mCi	X 1.00	= 107
Black & Red:	57.3 mCi	X 1.86	= 106.6
Black & Orange:	31.7 mCi	X 3.43	= 108.7
Black & Yellow:	7.8 mCi	X 13.67	= 106.6
Black & Green:	2.24 mCi	X 49.31	= 110.5
Black & Blue:	7.38 mCi	X 139.0	= 102.6
Black & Purple:	2.42 mCi	X 449.6	= 108.8
		SUM	= 750.8

$$\text{MEAN} = \frac{\text{SUM}}{7} = 107.3$$

$$\text{MEAN} \times 1.05 = 112.6 = \text{UPPER LIMIT}^*$$

$$\text{MEAN} \times 0.95 = 101.9 = \text{LOWER LIMIT}^*$$

Compare Column D data to upper and lower limits to confirm linearity.

*Instead of a variation in the Column D data of $\pm 6\%$, your radioactive material license may allow a difference of $\pm 10\%$ in the test results. If so, multipliers of 1.10 and 0.90 can be used to determine the upper and lower limits.

ACCURACY TEST

DOSE CALIBRATOR - DETROIT MEMORIAL - 5/31/85

SOURCE CALIB. DATE	A_0	CALCUL. A_N	DETECTED A_N	VARIANCE
CESIUM-137 8/4/78	227 uCi	194 uCi	192.7 uCi	-0.5%
BARIUM-133 8/4/78	207 uCi	126.4 uCi	121.0 uCi	-4.4%
COBALT-57	5.3 mCi	4.02 mCi	3.91 mCi	-2.9%

TESTS PERFORMED BY C. J. DAMICO, R.S.O.