



# ST. ELIZABETH HOSPITAL MEDICAL CENTER

A MAJOR TEACHING HOSPITAL ASSOCIATED WITH  
The Northeastern Ohio Universities College Of Medicine

1044 Belmont Avenue  
P.O. Box 1790  
Youngstown, Ohio 44501-1790

(216) 746-7211  
Extension 3666

DEPARTMENT OF  
RADIOLOGY

N. R. C. License # 34-01131-01

12 October, 1983

United States Nuclear Regulatory Commission  
Material Licensing Branch  
Region III  
799 Roosevelt Road  
Glenn Ellyn, Illinois 60137

Dear Sir:

Please amend our N. R. C. License # 34-01131-01 in three (3) areas as stated below.

1. Due to the increase in the volume of Nuclear Cardiology procedures, (3,000 procedures in 1982), please amend our current possession limit of Mo 99 generators from 2 Ci to 3 Ci.
2. We presently place our Mo99 generators in storage until radiation levels reach background. We wish to amend our license as to return spent generators to the manufacturer. Guide lines for this procedure will be followed according to N. R. C., Department of Transportation and Manufactures regulations.
3. Test for instrument linearity is presently being done on our dose calibrator in accordance with attachment No.#6., sub-stem E, of our original license application. (enclosure)

As an alternative to our present procedure, the dose calibrator can be checked for activity linearity with use of a device called "calicheck" form Calcorp, Inc. The manufacturers instructions for use as revised on March 2, 1982, will be followed. Test results will be recorded and retained for inspection. Corrective action will be taken if unacceptable linearity is demonstrated.

Thank you for your co-operation.

Sincerely,

*Clayton A. Hixson*

Clayton A. Hixson, M. D.

Medical Director--Division of Nuclear Medicine

Operated by the Sisters of the Humility of Mary

8506100018 850531

REG3 LIC30

34-01131-01

PDR

RECEIVED

OCT 17 1983

REGION III

030-02670

RECEIVED BY LFMB

Date: 10/21/83

Log: 10/21/83

By: [Signature]

Orig. To: [Signature]

Action Comm: [Signature]

11/2/83

Applicant: [Signature]

Check No. 119303

Amount/Fee Category \$407.40

Type of Fee: Amend

Date Check Rec'd: 11/18/83

Received By: [Signature]

OCT 17 1983

# CODE FOR CALIBRATION OF DOSE CALIBRATOR

## A. Test for the following:

1. Instrument linearity (at installation and after any major repair)
2. Geometrical variation (at installation)
3. Instrument accuracy (at installation and annually thereafter)

## B. After repair or adjustment of the dose calibrator, repeat all the appropriate tests listed above (dependent upon the nature of the repair).

## C. Daily or before each use of the instrument:

1. Measure and record the activity of at least one reference standard (e.g., 1-2 mCi of Co-57). This check should be repeated every day whenever sample readings are not within 10% of the expected assay. Variation greater than 5% in this test will indicate the need for instrument repair, adjustment, or recalibration.
2. Measure and record the apparent activity of a long-lived standard radionuclide such as Cs-137 or Ra-226. Choose a source with activity in the 100 uCi range.

## D. Inspect the instrument on a quarterly basis to ascertain that the instrument chamber liner is in place and that instrument zero is proper, etc. (see manufacturer's instructions).

## E. Test of Instrument Linearity

The linearity of a dose calibrator should be ascertained over the entire range of activities employed. This test will utilize a vial of Tc-99m whose activity is equivalent to the maximum anticipated activity to be assayed (e.g., the first elution from a new generator).

1. Assay the Tc-99m vial in the dose calibrator and subtract background level to obtain net activity in mCi/mCi.
2. Repeat step 1 at time intervals of 5, 20, 30, and 48 hours after the initial assay.

Using the 30-hour activity as a standard, the predicted activities at 0, 24, and 48 hours are as follows:

Assay Time (hr)	Correction Factor
0	12
6	14
24	3
30	1
48	0.125

Example: If the net activity measured at 30 hours was 1.025 mCi, the predicted activity for 0 and 48 hours would be 12.30 mCi and  $1.025 \text{ mCi} \times 0.125 = 0.128 \text{ mCi}$  respectively.

1. Plot the measured net activity for each time interval against the predicted activity on log-log graph paper.
2. The activities plotted should be within  $\pm 5\%$  of the predicted line if the instrument is linear and functioning properly. Errors greater than  $\pm 5\%$  indicate the need for repair or adjustment of the instrument.
3. If instrument linearity cannot be corrected, it will be necessary to perform assays to either assay an aliquot of the sample that was accurately measured or to use the graph constructed by the manufacturer to convert measured activities to true activities.

#### F. Test for Geometrical Variation

There may be significant geometrical variation in activity measured as a function of sample volume or configuration, depending on the volume and shape of the ionization chamber used in the assay. The extent of geometrical variation should be determined by comparing measured values with appropriate correction factors computed for variations in geometry. If, i.e., greater than  $\pm 2\%$  (or, if correction factors are not provided by the manufacturer, the accuracy of these should be  $\pm 2\%$ ).

To measure variation with volume of liquid, a 30-cc vial containing 10 cc of Co-57 or other appropriate radioactive isotope in a volume of 1 ml of liquid.

1. Assay vial at the appropriate instrument setting and correct for background level to obtain net activity.