

UNITED STATES ATOMIC ENERGY COMMISSION
APPLICATION FOR BYPRODUCT MATERIAL LICENSE

Form approved
Budget Bureau No. 38-80027

INSTRUCTIONS - Complete Items 1 through 16 if this is an initial application or an application for renewal of a license. Information contained in previous applications filed with the Commission with respect to Items 8 through 15 may be incorporated by reference provided references are clear and specific. Use supplemental sheets where necessary. Item 16 must be completed on all applications. Mail two copies to: U.S. Atomic Energy Commission, Washington, D.C. 20545, Attention: Isotopes Branch, Division of Materials Licensing. Upon approval of this application, the applicant will receive an AEC Byproduct Material License. An AEC Byproduct Material License is issued in accordance with the general requirements contained in Title 10, Code of Federal Regulations, Part 30, and the Licensee is subject to Title 10, Code of Federal Regulations, Part 20.

<p>1 (a) NAME AND STREET ADDRESS OF APPLICANT (Institution, firm, hospital, person, etc. Include ZIP Code)</p> <p>Department of the Army Fitzsimons Army Medical Center and U.S. Army Medical Research and Nutrition Laboratory Denver, Colorado 80240</p>	<p>(b) STREET ADDRESS(ES) AT WHICH BYPRODUCT MATERIAL WILL BE USED (If different from 1 (a) Include ZIP Code)</p> <p>Department of the Army Fitzsimons Army Medical Center and U.S. Army Medical Research and Nutrition Laboratory Denver, Colorado 80240</p>
<p>2 DEPARTMENT TO USE - BYPRODUCT MATERIAL</p> <p>Department of Radiology Nuclear Medicine Service</p>	<p>3 PREVIOUS LICENSE NUMBER(S) (If this is an application for renewal of a license, please indicate and give number)</p> <p>Amendment to Existing License No. 05-00046-13 (30 Apr 74)</p>
<p>4 INDIVIDUAL USER(S) (Name and title of individual(s) who will use or directly supervise use of byproduct material. Give training and experience in Items 8 and 9)</p> <p>See Application dated 12 March 1973 (Control No. 35871)</p>	<p>5 RADIATION PROTECTION OFFICER (Name of person designated as radiation protection officer if other than individual user. Attach resume of his training and experience as in Items 8 and 9)</p> <p>See Application dated 12 March 1973 (Control No. 35871)</p>
<p>6 (a) BYPRODUCT MATERIAL: (Elements and mass number of each)</p> <p>Molybdenum-99</p>	<p>(b) CHEMICAL AND OR PHYSICAL FORM AND MAXIMUM NUMBER OF MILLCURIES OF EACH CHEMICAL AND OR PHYSICAL FORM THAT YOU WILL POSSESS AT ANY ONE TIME (If sealed source(s), also state name of manufacturer, model number, number of sources and maximum activity per source)</p> <p>Contained in the E.R. Squibb and Sons, Inc., Minitec sterile Tc-99m generator as a replacement for the generator E.R. Squibb and Sons, Inc. Model 08871 which will no longer be manufactured as of 1 May 1974.</p> <p>1 curie as currently authorized.</p> <p>or - in any AEC approved Tc-99m sterile generator recommended by the AEC as a substitute.</p>
<p>7 DESCRIBE PURPOSE FOR WHICH BYPRODUCT MATERIAL WILL BE USED (If byproduct material is for human use, Supplement A (Form AEC 313a) must be completed in lieu of this item. If byproduct material is in the form of a sealed source, include the make and model number of the storage container and or device in which the source will be stored and or used.)</p> <p>The generator is used as a source of Tc-99m pertechnetate for subsequent use in a number of standard diagnostic procedures authorized under the license.</p> <p>The generator will be used in accordance with manufacturer's instructions. Nuclear Medicine Personnel will be instructed in the health and safety precautions involved in its use by the Chief of Nuclear Medicine and the Radiation Protection Officer.</p> <p>See the attached manufacturer's instruction summary sheet.</p>	

46485

TRAINING AND EXPERIENCE OF EACH INDIVIDUAL NAMED IN ITEM 4 (Use supplemental sheets if necessary)

B. TYPE OF TRAINING	WHERE TRAINED	DURATION OF TRAINING	ON THE JOB (Circle answer)	FORMAL COURSE (Circle answer)
a. Principles and practices of radiation protection	See application dated 12 March 1974 (Control No. 35871)		Yes No	Yes No
b. Radioactivity measurement standardization and monitoring techniques and instruments			Yes No	Yes No
c. Mathematics and calculations basic to the use and measurement of radioactivity			Yes No	Yes No
d. Biological effects of radiation			Yes No	Yes No

9. EXPERIENCE WITH RADIATION (Actual use of radionuclides or equivalent experience)

ISOTOPE	MAXIMUM AMOUNT	WHERE EXPERIENCE WAS GAINED	DURATION OF EXPERIENCE	TYPE OF USE
See application dated 12 March 1974 (Control No. 35871)				

10. RADIATION DETECTION INSTRUMENTS (Use supplemental sheets if necessary)

TYPE OF INSTRUMENTS (Include name and model number of each)	NUMBER AVAILABLE	RADIATION DETECTED	SENSITIVITY RANGE (mr/hr)	WINDOW THICKNESS (mg/cm)	USE (Monitoring, surveying, measuring)
See application dated 12 March 1973 (Control No. 35871) and application dated 25 June 1968.					

11. METHOD, FREQUENCY, AND STANDARDS USED IN CALIBRATING INSTRUMENTS LISTED ABOVE

See application dated 12 March 1973 (Control No. 35871) and application dated 25 June 1968.

12. FILM BADGES, DOSIMETERS, AND BIO ASSAY PROCEDURES USED (For film badges, specify method of calibrating and processing, or name of supplier)

See application dated 12 March 1973 (Control No. 35871) and application dated 25 June 1968.

INFORMATION TO BE SUBMITTED ON ADDITIONAL SHEETS IN DUPLICATE

13. FACILITIES AND EQUIPMENT Describe laboratory facilities and equipment handling equipment, storage containers, shielding, work habits, etc. (Sketch of facility if checked) (Circle answer) Yes No
See application dated 12 March 1973 (Control No. 35871) and application dated 25 June 1968.

14. RADIATION PROTECTION PROGRAM Describe the radiation protection program including control measures. If application covers sealed sources, submit leak test records, if applicable, for all sealed sources. If application covers unsealed sources, submit records of contamination monitoring, maintenance and repair of the source.
See application dated 12 March 1973 (Control No. 35871) and application dated 25 June 1968.

15. WASTE DISPOSAL Describe the disposal of radioactive waste. Specify the type and amount of waste, the method of disposal, and the location of the disposal site. If application covers sealed sources, submit records of leak testing and disposal of the source.
See application dated 12 March 1973 (Control No. 35871) and application dated 25 June 1968.

CERTIFICATE (This item must be completed by applicant)

16. THE APPLICANT AND ANY OFFICIAL EXECUTING THIS CERTIFICATE ON BEHALF OF THE APPLICANT NAMED IN ITEM 1, CERTIFY THAT THIS APPLICATION IS PREPARED IN CONFORMITY WITH TITLE 10, CODE OF FEDERAL REGULATIONS, PART 30, AND THAT ALL INFORMATION CONTAINED HEREIN, INCLUDING ANY SUPPLEMENTS ATTACHED HERETO, IS TRUE AND CORRECT TO THE BEST OF OUR KNOWLEDGE AND BELIEF.

Fittsimons Army Medical Center
Denver, Colorado 80240

By H. F. Cowgill, M.D.
Colonel, MC

Date _____

Title of certifying official

WARNING.—18 U. S. C., Section 1001, Act of June 25, 1948, 62 Stat. 749, makes it a criminal offense to make a willfully false statement or representation to any department or agency of the United States as to any matter within its jurisdiction.



CAUTION: NEW DRUG—Limited by United States law to investigational use

MINITECTM

Technetium 99m

GENERATOR

Minitec (Technetium 99m) Generator provides a means of obtaining a sterile, non-pyrogenic supply of Technetium 99m (^{99m}Tc), as sodium (^{99m}Tc) pertechnetate, a versatile scanning agent that can be administered intravenously or orally. ^{99m}Tc, the short-lived daughter ($T_{1/2}=6$ hours) of Molybdenum 99 (⁹⁹Mo, $T_{1/2}=67$ hours), is obtained from the generator by periodic elution. The amount (in millicuries) of ^{99m}Tc obtained in the initial elution will depend on the original potency of the generator, while the activity obtained from subsequent elutions will depend on the time interval between elutions.

Eluting the generator every 24 hours will provide optimal amounts of ^{99m}Tc. Most laboratories will therefore find it convenient to elute the generator each day at a specific time. However, the generator may be eluted whenever sufficient amounts of ^{99m}Tc have accumulated within the column.

DESCRIPTION

The fission product Molybdenum 99 used in the generator meets or exceeds the purity requirements of the Atomic Energy Commission with respect to allowable levels of contaminants. The aluminum concentration is not more than 10 mcg. per cc. of generator eluate. The eluate meets the A.E.C. limits for ⁹⁹Mo contamination, that is, not more than one microcurie ⁹⁹Mo per millicurie of ^{99m}Tc or five microcuries per dose of ^{99m}Tc administered. Each elution from the generator should be assayed before use for ^{99m}Tc activity and for the possible presence of ⁹⁹Mo. Directions for both assays are provided in this monograph. Material containing more than 5 microcuries of ⁹⁹Mo per dose of ^{99m}Tc pertechnetate exceeds Atomic Energy Commission limits and should not be administered. The generator consists of a specially designed lead shield containing an alumina-packed column which releases ^{99m}Tc upon elution. The lead shield has access ports at the top of the column, allowing aseptic elution and storage under conditions of constant shielding.

Supplied with the generator are vials of sterile, non-pyrogenic eluent, and suitable equipment for eluting, collecting, and assaying the Technetium 99m.

WARNINGS

Due to the high precalibration activities of technetium generators, the individual user should make certain that he takes proper precautions to insure that he is within his permitted possession limits (e.g. a 300 mCi generator calibrated for noon Friday will have an activity of approximately 845 mCi at 8 A.M. on the preceding Monday). In order to evaluate radiation hazards properly in the event of accidents or spills, it is essential that personnel be aware of the actual quantity of activity involved.

Maintain proper radiation safety precautions at all times. The column containing ⁹⁹Mo must not be removed from the lead shield at any time. The radiation field surrounding an unshielded column is quite high. Solutions of ^{99m}Tc withdrawn from the generator should always be adequately shielded. The early elutions from the generator are highly radioactive.

IMPORTANT

Since material obtained from the generator may be intended for intravenous administration, aseptic technique must be strictly observed in all handling. Only the eluent provided should be used to elute the generator. Do not administer material eluted from the generator if there is any evidence of foreign matter.

DIRECTIONS FOR ELUTING ^{99m}Tc

Read entire procedure before beginning elution. Carefully follow each step in the order described.


1. Attach sterile Needle Adapter Assembly to the lead generator shield.
 - (a) Remove the two port plugs from top of the shield.
 - (b) Using a sterile cotton-tipped applicator, swab the exposed rubber closures inside the port holes with a suitable germicide.
 - (c) Position Assembly so that the guide pin on underside of Assembly will align with hole in shield; the needle ends will then align with port holes. (Assembly can only be mounted one way.)
 - (d) Press Assembly into port holes so that it is firmly seated. The needles have now pierced the rubber closures. Be certain that Assembly is securely attached to the shield. Once seated, the Assembly should not be removed.
2. Attach Sterile Eluent vial:
(Fits into the short compartment on Needle Adapter Assembly)
 - (a) Swab rubber closure of eluent vial with germicide.
 - (b) Remove cap from the short compartment on Assembly and firmly position eluent vial on the needles.

NOTE: THE ELUENT VIAL MUST BE ATTACHED TO THE NEEDLE ADAPTER ASSEMBLY BEFORE ATTACHING THE EVACUATED VIAL. OR THE VACUUM IN THE EVACUATED VIAL WILL BE LOST.

48405 J3-281

3. Attach Sterile Evacuated Collecting Vial:
(Fits into the tall compartment on Needle Adapter Assembly)
 - (a) Place collecting vial into the Collecting Vial lead shield; cover with the lead neck retainer cover and firmly secure cover to shield with the plastic lock ring.
 - (b) Swab rubber closure of the lead-shielded collecting vial with germicide
 - (c) Remove cap from the tall compartment on Assembly and firmly position collecting vial on the needle. This will start the elution. Elution should be complete 3 minutes after eluent vial has emptied.
4. When elution is complete, replace collecting vial with the Minitec Needle Cover vial.
 - (a) Swab rubber closure of needle cover vial with germicide.
 - (b) Remove lead-shielded collecting vial (which now contains the radioactive eluate) from the needle and immediately position needle cover vial on the needle. Push needle cover vial on needle as far as it will go.
 - (c) For added shielding of eluate, place lead cap over the lead-shielded collecting vial.
5. To maintain sterility of the system, after each milking, leave the empty Sterile Eluent vial (Step 2) in position and replace the collecting vial with the same needle cover vial (Step 4) until the generator is to be eluted again.

The following illustration shows generator assembled for elution.



DIRECTIONS FOR ASSAYING ^{99m}Tc ACTIVITY

The ^{99m}Tc activity of the Minitec Generator eluate may be assayed by the Whole Vial Assay Method described below. (Note: The Minitec Whole Vial Assay Kit is available upon request only.)

A ^{99m}Tc standard contained in a 5 cc. volume should be used to establish the ^{99m}Tc factor. A container having the same geometry as the collecting vial should be used to contain the 5 cc. ^{99m}Tc standard.

The Whole Vial Assay method can be performed with any scintillation detector used in conjunction with a pulse height analyzer (well-type scintillation spectrometer or an upright scintillation probe). The Whole Vial Assay components are: (A) top lid, (B) ring lid, (C) shield, (D) Cobalt 57 check source (from 5 to 15 μCi), (E) platform, and (F) lead absorbers.

1. Place the shield (C) on top of the well of a scintillation spectrometer; if the well aperture is too large for support, place the shield in the middle of the platform (E) and center over the well.
2. Set the discriminator to the 100 Kev to 200 Kev range. Record the background counting rate.
3. Insert a ^{99m}Tc standard* in the shield (C) and cover it with ring lid (B) and the top lid (A).
4. Count and record the ^{99m}Tc standard activity. (NOTE: If the count rate is too high, place lead absorbers (F) into the shield (C) until an acceptable count rate is obtained. To prevent errors, the lead absorbers should then be glued to the bottom of the shield.)
5. Replace the ^{99m}Tc standard with the ^{99m}Tc sample, cover as in step 3, and record the count rate of the ^{99m}Tc sample. (NOTE: The volume of the standard solution and the sample solution must be equal and in a container of the same geometry.)
6. Replace the shield (C) with the Cobalt 57 check source (D) and record count rate. [NOTE: Once the count rate of the Cobalt 57 check source has been determined, the continued use of the ^{99m}Tc standard is unnecessary provided the count rate of the check source does not vary more than $\pm 10\%$ from the expected count rate (after correction for radioactive decay).]
7. For subsequent assays only steps 5 and 6 need be followed.
8. Calculate the ^{99m}Tc factor as follows:

$$^{99m}\text{Tc factor} = \frac{\text{activity (mCi) of } ^{99m}\text{Tc standard}}{\text{net count rate (cpm) of } ^{99m}\text{Tc standard}}$$

Example:

activity of ^{99m}Tc standard = 10 mCi
count rate as determined in step 4 = 20,000 cpm

$$^{99m}\text{Tc factor} = \frac{10 \text{ mCi}}{20,000 \text{ cpm}} = 0.0005$$

9. The ^{99m}Tc factor (step 8) is used to determine the ^{99m}Tc activity in 5 cc. of eluate from the Sterile Generator provided the count rate of the Cobalt 57 check source is within $\pm 10\%$ of the expected count rate. Calculate as follows:

$$^{99m}\text{Tc activity (mCi) of } ^{99m}\text{Tc sample} = ^{99m}\text{Tc factor} \times \text{net count rate (cpm) of } ^{99m}\text{Tc sample}$$

*A ^{99m}Tc sample calibrated by existing methods, such as the dilution method, can be used as the initial standard.

Example:

$$^{99m}\text{Tc factor} = 0.0005$$

$$\text{net count rate from unknown sample} = 38,000 \text{ cpm}$$

$$^{99m}\text{Tc activity (mCi/cc)} = 0.0005 \times 38,000 = 19 \text{ mCi/cc}$$

TEST FOR INSTRUMENT LINEARITY

The use of ^{99m}Tc generators has resulted in the assay of considerably larger quantities of activity (up to 1 Ci of ^{99m}Tc) than previously encountered. It would appear worthwhile, therefore, to check the accuracy of the assay methods for these larger quantities of activity. This is especially important since most instruments are calibrated with standards containing approximately 10 mCi of ^{99m}Tc while the initial eluates contain up to 900 mCi of ^{99m}Tc . It is not unlikely that spectrometers used in the whole vial assay procedure are not linear up to the count rates being encountered in initial assays, and that dose calibrators standardized on the lower scale are not accurate on the higher scales.

To check the assay procedures where yield problems are encountered with the high activity generators, the following is recommended:

1. Assay the entire ^{99m}Tc eluate on the first day of generator use by the procedure normally used. Record assay in millicuries/cc.
2. Withdraw 1 cc. of the eluate and transfer to an empty collecting vial. Add 4 cc. of water and mix thoroughly.
3. Assay the collecting vial containing the eluate-water solution. Record assay, corrected for decay, in total millicuries present in the vial.
4. If the calibration of the assay procedure is linear for the higher activities, the values obtained in Steps 1 and 3 should be identical. If the instrument is losing counts at the higher activities, the value obtained in Step 3 will be higher than that obtained in Step 1.

For example:

- (a) The assay value in Step 1 is 660 mCi in 5 cc., therefore, the eluate contains 132 mCi/cc.
- (b) The assay value in Step 3 is 169 mCi.
- (c) Assuming that the 169 mCi value is correct, since the instrument was calibrated in this range of activity, the assay would be off (on the low side) by approximately 22%. Calculated as follows:

$$\begin{aligned}\% \text{ error} &= \frac{\text{True activity} - \text{Observed activity}}{\text{True activity}} \times 100 \\ &= \frac{169 - 132}{169} \times 100 \\ &= 22\%\end{aligned}$$

If an error of greater than 5% is observed in conducting the above assay check, it will be necessary in routine assays to either assay an aliquot of the eluate that can be accurately measured, or construct a curve to correct for counting losses at higher counting rates. Such a curve can be constructed by assaying various aliquots of a high activity milking—e.g. 1, 2, 3, 4, and 5 cc.—and plotting the activity assayed versus the true activity (calculated for assay of 1 cc. aliquot).

DIRECTIONS FOR ASSAYING ^{99}Mo ACTIVITY

A. Whole Vial Assay Method

1. Determination of the ^{99}Mo factor:

Use the Whole Vial Assay components previously described and set the discriminator to the 0.6 Mev to 1.0 Mev range. Calculate the ^{99}Mo factor by dividing the activity (in μCi) of the Cesium 137 Standard by the net count rate (gross cpm minus background cpm) of the Cesium 137 Standard and multiplying by a factor of 4.5 to convert ^{137}Cs activity to equivalent ^{99}Mo activity. (NOTE: The Cesium 137 Standard is provided on request in a collecting vial with a total volume of 5 cc. If the Squibb Standard is not used, be certain that the substitute ^{137}Cs standard is used in a container having the same geometry as the collecting vial.)

Calculate as follows:

$$^{99}\text{Mo factor} = \frac{\text{Cesium 137 activity (total } \mu\text{Ci}) \times 4.5}{\text{net count rate of Cesium 137 Standard (cpm)}}$$

NOTE: The ^{99}Mo factor remains the same for each assay if the instrument calibration is not altered; if there is a change in calibration, a new factor must be determined.

2. For each assay:

- a. Place the lead shield on top of the well of a scintillation spectrometer; if the well aperture is too large for support, place the shield in the middle of the Whole Vial Assay platform and center over the well.
- b. Set the discriminator to the 0.6 Mev to 1.0 Mev range. Record the background activity (cpm).
- c. Insert the vial containing the ^{99m}Tc sample into the lead shield and cover with the ring and top lids.
- d. Record the net count rate (cpm) of the ^{99m}Tc sample (gross cpm minus background cpm).
- e. Calculate the ^{99}Mo activity as follows:

$$^{99}\text{Mo activity} = ^{99}\text{Mo factor (from step 1)} \times \text{net count rate (cpm) of } ^{99m}\text{Tc sample}$$

Example:

$$\begin{aligned}\text{net count rate of } ^{99m}\text{Tc sample} \\ (\text{as determined in step 2 d}) &= 1700 \text{ cpm} \\ ^{99}\text{Mo factor} &= 0.003 \mu\text{Ci/cpm} \\ ^{99}\text{Mo activity} &= 0.003 \mu\text{Ci/cpm} \times 1700 \text{ cpm} \\ &= 5.1 \mu\text{Ci}\end{aligned}$$

- f. Calculate the concentration of ^{99}Mo per millicurie of ^{99m}Tc sample by dividing the ^{99}Mo activity (from step 2 e) by the total ^{99m}Tc activity.

Example:

$$\frac{^{99}\text{Mo activity (in } \mu\text{Ci})}{^{99m}\text{Tc activity (in mCi)}} = \frac{5.1}{51} = 0.1 \mu\text{Ci Molybdenum } 99 \text{ per mCi Technetium } 99m$$

B. Alternate Method

1. Place the collecting vial containing the total ^{99m}Tc eluate in a $\frac{1}{4}$ -inch lead container and set this on the surface of a well-type scintillation detector or scintillation probe. (The lead container is provided on request with the Cesium 137 Standard mentioned below.)
2. Determine the activity in the 0.6—1.0 Mev range. Record net counts/minute.

