



March 2, 1984

Procedure #26Vacuuuming of External Block of Old Cell

Vacuuuming of external block of old cell will be carried out by the following procedure:

- A.) Safety of Personnel: Same as Procedure #23.
- B.) Work to be Performed:
 - 1.) All external block will be vacuuumed. All material picked up by vacuum will be placed into a plastic bag. This material will be surveyed & placed into a drum designated DSM-55 LSA.
 - 2.) Vacuum will be wiped down & will be left inside enclosure to be used in later procedures.
- C.) Tools:
 - 1.) Sears wet/dry electric vacuum cleaner
 - 2.) Eberline 130-A

JJD/lm

8506180269 850426
PDR FOIA
TERPILAB4-763 PDR

296



March 2, 1984

Procedure #26

Vacuuming of External Block of Old Cell

Vacuuming of external block of old cell will be carried out by the following procedure:

- A.) Safety of Personnel: Same as Procedure #23.
- B.) Work to be Performed:
 - 1.) All external block will be vacuumed. All material picked up by vacuum will be placed into a plastic bag. This material will be surveyed & placed into a drum designated DSM-55 LSA.
 - 2.) Vacuum will be wiped down & will be left inside enclosure to be used in later procedures.
- C.) Tools:
 - 1.) Sears wet/dry electric vacuum cleaner
 - 2.) Eberline 130-A

JJD/lnm